

Product Manual

96- Well Ras Activation ELISA Kit (Chemiluminescent)

Catalog Number

A6601, 96T

96 assays

MyBioSource.com

Introduction

Small GTP-binding proteins (or GTPases) are a family of proteins that serve as molecular regulators in signaling transduction pathways. Ras, a 21 kDa protein, regulates a variety of biological response pathways that include cell growth, cell transformation and tumor invasion. Like other small GTPases, Ras regulates molecular events by cycling between an inactive GDP-bound form and an active GTP-bound form. In its active (GTP-bound) state, Ras binds specifically to the Ras-binding domain (RBD) of Raf-1 to control downstream signaling cascades. The most notable members of the Ras subfamily are H-Ras, N-Ras and K-Ras, mainly for being implicated in many types of cancer.

96-well Ras Activation ELISA Kit utilizes plate-bound, Raf-1 RBD to selectively isolate and pull-down the active forms of Ras (H-, K-, and N-Ras isoforms from human, mouse and rat) from purified samples or endogenous lysates. Subsequently, the captured GTP-Ras is detected by an Anti-pan-Ras Antibody and HRP conjugated secondary antibody. Finally, following the addition of Chemiluminescent Reagent, the reaction is then measured in a plate luminometer.

96-well Ras Activation ELISA Kit provides a simple and fast tool to monitor the activation of Ras. Each kit provides sufficient reagents to perform up to 96 assays.

Related Products

1. H-Ras Activation Assay Kit
2. K-Ras Activation Assay Kit
3. N-Ras Activation Assay Kit
4. Rac1 Activation Assay Kit
5. Cdc42 Activation Assay Kit
6. RhoA Activation Assay Kit
7. Raf-1 RBD Agarose Beads
8. 96-well Ras Activation ELISA Kit (Colorimetric)
9. Ras Expression Vector Set
10. Active Ras Expression Vector Set

Kit Components

Box 1 (shipped at room temperature)

1. Raf-1 RBD Capture Plate (Part No. 244101): One solid, white 96-well plate.
2. 5X Assay/Lysis Buffer (Part No. 240102): One 30 mL bottle of 125 mM HEPES, pH 7.5, 750 mM NaCl, 5% NP-40, 50 mM MgCl₂, 5 mM EDTA, 10% Glycerol.
3. Assay Diluent (Part No. 310804): One 50 mL bottle.
4. 10X Wash Buffer (Part No. 310806): One 100 mL bottle.
5. Chemiluminescent Reagent A (Part No. 250102): One 6 mL amber bottle.

6. Chemiluminescent Reagent B (Part No. 250103): One 6 mL bottle.
7. Plate Sealing Films (Part No. 244102): One pack of two sealing films.

Box 2 (shipped on blue ice packs)

1. Raf-1 RBD (500X) (Part No. 244002): One 40 μ L vial.
2. Anti-pan-Ras Antibody (1000X) (Part No. 244003): One 20 μ L vial.
3. Secondary Antibody, HRP Conjugate (Part No. 244004): One 20 μ L vial.
4. 100X GTP γ S (Part No. 240103): One 50 μ L vial of 10 mM GTP γ S dissolved in sterile water.
5. 100X GDP (Part No. 240104): One 50 μ L vial of 100 mM GDP dissolved in sterile water.

Materials Not Supplied

1. Stimulated and non-stimulated cell or tissue lysates
2. Ras activators
3. Protease inhibitors
4. 0.5 M EDTA in water
5. 1 M MgCl₂
6. 30°C incubator or water bath
7. Room temperature shaker
8. 10 μ L to 1000 μ L adjustable single channel micropipettes with disposable tips
9. 50 μ L to 300 μ L adjustable multichannel micropipette with disposable tips
10. Multichannel micropipette reservoir
11. 96-well Plate Luminometer

Storage

Upon receipt, aliquot and store Raf-1 RBD at -80°C and avoid freeze/thaw. Aliquot and store the Anti-pan-Ras Antibody, GTP γ S, and GDP components at -20°C and avoid freeze/thaw. Store all other components at 4°C.

Preparation of Reagents

- 1X Assay/Lysis Buffer: Mix the 5X stock briefly and dilute to 1X in deionized water. Just prior to usage, add protease inhibitors such as 1 mM PMSF, 10 μ g/mL leupeptin, and 10 μ g/mL aprotinin.
- 1X Wash Buffer: Dilute the 10X Wash Buffer to 1X with deionized water. Stir to homogeneity.
- Anti-pan-Ras Antibody: Immediately before use dilute the Anti-pan-Ras Antibody 1:1000 with Assay Diluent. Do not store diluted solutions.
- Secondary Antibody, HRP Conjugate: Immediately before use dilute the Secondary Antibody, HRP Conjugate 1:2500 with Assay Diluent. Do not store diluted solutions.

- Chemiluminescent Reagent: Immediately before use, mix equal volumes of Chemiluminescent Reagent A with Chemiluminescent Reagent B. Do not store diluted solutions.

Preparation of Samples

Note: It is advisable to use fresh cell or tissue lysates because GTP-Ras is quickly hydrolyzed to GDP-Ras; frozen lysates stored at -70°C may be used. Performing steps at 4°C or on ice may reduce hydrolysis. Avoid multiple freeze/thaw cycles of lysates.

I. Adherent Cells

1. Culture cells to approximately 80-90% confluence. Stimulate cells with Ras activator(s) as desired.
2. Aspirate the culture media and wash twice with ice-cold PBS.
3. Completely remove the final PBS wash and add ice-cold 1X Assay/Lysis Buffer to the cells (0.5 - 1 mL per 100 mm tissue culture plate).
4. Place the culture plates on ice for 10-20 minutes.
5. Detach the cells from the plates by scraping with a cell scraper.
6. Transfer the lysates to appropriate size tubes and place on ice.
7. If nuclear lysis occurs, the cell lysates may become very viscous and difficult to pipette. If this occurs, lysates can be passed through a 27½-gauge syringe needle 3-4 times to shear the genomic DNA.
8. Clear the lysates by centrifugation for 10 minutes (14,000 x g at 4°C).
9. Collect the supernatant and store samples on ice for immediate use, or snap freeze samples and store at -70°C for future use.
10. Proceed to GTPγS/GDP Loading for positive and negative controls, or the Activation ELISA (Assay Protocol Section).

II. Suspension Cells

1. Culture cells and stimulate with Ras activator(s) as desired.
2. Perform a cell count, and then pellet the cells by centrifugation.
3. Aspirate the culture media and wash twice with ice-cold PBS.
4. Completely remove the final PBS wash and add ice-cold 1X Assay/Lysis Buffer to the cell pellet (0.5 – 1 mL per 1×10^7 cells).
5. Lyse the cells by repeated pipetting.
6. Transfer the lysates to appropriate size tubes and place on ice.
7. If nuclear lysis occurs, the cell lysates may become very viscous and difficult to pipette. If this occurs, lysates can be passed through a 27½-gauge syringe needle 3-4 times to shear the genomic DNA.
8. Clear the lysates by centrifugation for 10 minutes (14,000 x g at 4°C).

9. Collect the supernatant and store samples on ice for immediate use, or snap freeze samples and store at -70°C for future use.
10. Proceed to GTP γ S/GDP Loading for positive and negative controls, or the Activation ELISA (Assay Protocol Section).

Assay Protocol

I. GTP γ S/GDP Loading (Positive and Negative Controls)

Note: Samples that will not be GTP γ S/GDP loaded may be kept on ice during preparation of GTP γ S/GDP loading samples.

1. Aliquot 0.5 mL of each cell lysate to two microcentrifuge tubes.
Note: Typical protein concentration of sample is > 0.5 mg/mL.
2. Add 10 μ L of 0.5 M EDTA to each sample.
3. Add 5 μ L of 100X GTP γ S to one tube (positive control) and 5 μ L of 100X GDP to the other tube (negative control). Mix and label each tube appropriately.
4. Incubate the tubes for 30 minutes at 30°C with agitation.
5. Stop the loading by adding 33 μ L of 1 M MgCl₂ to each tube. Mix and place tubes on ice.
6. Continue with the Activation ELISA.

II. Ras Activation ELISA

Note: Samples and controls should be thawed/maintained on ice just prior to use (Step 3).

1. Determine the number of wells to be used, and dilute the Raf-1 RBD 1:500 in Assay Diluent. Add 100 μ L of the diluted Raf-1 RBD to each well of the Raf-1 RBD Capture Plate. Incubate at room temperature for 1 hour on an orbital shaker. Cover unused wells with a plate sealing film.
Note: Do not store diluted solutions.
2. Wash microwell strips 3 times with 250 μ L 1X Wash Buffer per well with thorough aspiration between each wash. After the last wash, empty wells and tap microwell strips on absorbent pad or paper towel to remove excess 1X Wash Buffer.
3. Add 50 μ L of Ras lysate sample (10-100 μ g), control, or buffer blank per well. Each sample should be assayed in duplicate. Any sample dilutions should be performed in cold, 1X Assay/Lysis Buffer.
4. Immediately add 50 μ L of Assay Diluent to each well (100 μ L total volume). Incubate at room temperature for 1 hour on an orbital shaker.

5. Wash microwell strips 5 times with 250 μ L 1X Wash Buffer per well with thorough aspiration between each wash. After the last wash, empty wells and tap microwell strips on absorbent pad or paper towel to remove excess 1X Wash Buffer.
6. Add 100 μ L of the diluted Anti-pan-Ras Antibody to each well. Incubate at room temperature for 1 hour on an orbital shaker.
7. Wash the strip wells 5 times according to step 5 above.
8. Add 100 μ L of the diluted Secondary Antibody, HRP Conjugate to each well. Incubate at room temperature for 1 hour on an orbital shaker.
9. Wash the strip wells 5 times according to step 5 above. Proceed immediately to the next step.
10. Add 100 μ L of Chemiluminescent Reagent (see Preparation of Reagents Section) to each well, including the blank wells. Incubate at room temperature for 5 minutes on an orbital shaker.
11. Read the luminescence of each microwell on a plate luminometer.

Example of Results

The following figure demonstrates typical results seen with 96-well Ras Activation ELISA Kit. One should use the data below for reference only.




Figure 1: EGF Stimulation. HeLa cells were serum starved for 18 hours before EGF stimulation (50 ng/mL for 2 minutes). Lysates were then prepared according to Assay Protocol. Background has been subtracted from data.




Figure 2: Pan-Ras Antibody Specificity. Anti-pan-Ras Antibody specificity to purified, H-, K-, and N-Ras human isoforms by dot blot.

References

1. Bar-Sagi D., and Hall A. (2000) *Cell* **103**: 227-38.
2. de Rooij J., and Bos J. L. (1997) *Oncogene* **14**: 623-5.

Recent Product Citation

Kim, S. M. et al. (2016). EGFR-mediated reactivation of MAPK signaling induces acquired resistance to GSK2118436 in BRAF V600E mutant NSCLC cell lines. *Mol Cancer Ther.* doi:10.1158/1535-7163.MCT-15-0375.

MyBioSource.com